



ELISA 20× Washing Buffer

Product information

Catalog Number: CEB-C301

Size: 50 mL/bottle

ELISA 20× Washing Buffer are solutions used to eliminate unnecessary and unbound components that can cause background detection and interference in an ELISA assay. The buffer is the same formulation that is supplied with most ClinMax™ ELISA kits for cytokines and biomarkers, but not all. Please contact Technical Assistance for kit-specific information.

Storage

Keep the reagents stored at 2-30 °C. For opened kit and reconstituted reagents, can be stored for up to 30 days at 2-8 °C.

Prepare the working buffers.

1×Washing Buffer: Dilute 50 mL ELISA 20×Washing Buffer with ultra purified or deionized or distilled water to 1000 mL.

NOTE:

This kit is for research use only and is not for use in diagnostic or therapeutic applications.

Please do not use the reagents after the expiration date indicated on the label.

Manufactured and distributed by
ACRODiagnostics Inc.
TEL: 010-67855298-8147

Distributed by:
ACROBiosystems Inc.
US & Canada TEL: +1 800-810-0816
Asia & Pacific TEL: +86 400-682-252

ELISA 20X Washing Buffer User guide