

Mouse ROR2 / NTRKR2 Protein, His Tag (MALS verified)

Catalog # RO2-M52H3



BIOSYSTEMS
Acro
Surprise Inside!

Synonym

ROR2,NTRKR2

Source

Mouse ROR2 Protein, His Tag(RO2-M52H3) is expressed from human 293 cells (HEK293). It contains AA Glu 34 - Gly 403 (Accession # [Q9Z138-1](#)).

Predicted N-terminus: Glu 34

Molecular Characterization

ROR2(Glu 34 - Gly 403)
Q9Z138-1 Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 43.1 kDa. The protein migrates as 55-60 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Mouse ROR2 Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS



The purity of Mouse ROR2 Protein, His Tag (Cat. No. RO2-M52H3) is more than 85% and the molecular weight of this protein is around 45-65 kDa verified by SEC-MALS.

[Report](#)

Background

Tyrosine-protein kinase transmembrane receptor ROR2 is also known as Neurotrophic tyrosine kinase, receptor-related 2 (NTRKR2), which belongs to the protein kinase superfamily and Tyr protein kinase family and ROR subfamily. ROR2 is a homodimer protein, which can binds YWHAB or interact with WTIP. ROR2 may be involved in the early formation of the chondrocytes. It seems to be required for cartilage and growth plate development.

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

7/18/2024



Clinical and Translational Updates

**Discounts, Gifts,
and more!**

